Wayen Exosome Isolation Kit

User Manual


Cat# EIQ3-01001 (Serum) Version 2017-01

Description (For research only)

EIQ3-01001 kit is used to isolate / purify exosomes between 30 and 200 nm diameter from serum. By adding appropriate amount of reagents (A and B) to serum sample, and incubating within a short period, the exosomes can be collected after centrifugation.

Advantages

- ✓ Quantity: Higher yield (versus other kits or methods)
- ✓ Quality: Purer exosome (less serum high-abundant protein)
- ✓ Quick: Faster (< 2 hours)

Contents

EIQ3-01001 kit contains Reagent A and Reagent B.

Storage

The kits are shipped at room temperature and should be stored at 4 °C after received. Properly stored kits are stable for 1 years from the date received.

Components	Storage	Amount
Reagent A	4 °C	7.5 mL
Reagent B	4 °C	7.5 mL

Experiment Protocol of EIQ3-01001 (Serum)

1. Prepare Serum Sample

1.1 Take the serum sample from storage and keep it on ice. If starting with frozen sample, thaw the sample completely in a 25 °C water bath and then place it on ice.

- 1.2 Centrifuge the serum sample at $3000 \times g$ for 10 minutes at $4 \text{ }^{\circ}\text{C}$.
- 1.3 Transfer the supernatant to a new tube and place it on ice until ready to perform the isolation.

2. Isolate Exosomes (Balance the Reagent A and Reagent B to room temperature before use and the starting volume of serum is recommended to be 200 μL . The example below is shown with 200 μL serum)

- 2.1 Take out 200 μL pre-treated serum. Add 50 μL Reagent A to the serum and mix well by pipetting up and down or vortexing until obtain a homogenous mixture.
※ **Note: Serum : Reagent A = 4:1 (volume)**
- 2.2 Incubate the mixture at $4 \text{ }^{\circ}\text{C}$ for 30 minutes.
※ **Note: Do not rotate or shake the tube during the incubation period**
- 2.3 After incubation, centrifuge at $3000 \times g$ for 10 minutes at room temperature. Remove the supernatant and the pellet is at the bottom of the tube.
- 2.4 Centrifuge the pellet within a short moment and remove the residual supernatant.
- 2.5 Resuspend the pellet completely by adding 200 μL $1 \times$ sterile PBS. Mix well by pipetting up and down or vortexing until obtain a homogenous mixture.
※ **Note: Volume of re-suspension is equal to the starting serum volume at this procedure**
- 2.6 Add 50 μL Reagent B to the re-suspension and mix well by pipetting up and down or vortexing until obtain a homogenous solution.
※ **Note: The re-suspension : Reagent B = 4:1 (volume)**
- 2.7 Incubate the mixture at $4 \text{ }^{\circ}\text{C}$ for 30 minutes.
- 2.8 After incubation, centrifuge the mixture at $3000 \times g$ for 10 minutes at room temperature and remove the supernatant.
- 2.9 Centrifuge the sample within a short moment and remove the residual supernatant.
- 2.10 Resuspend the exosome pellet completely in 50-120 μL $1 \times$ sterile PBS and mix well to obtain a homogenous mixture. Once the pellet is re-suspended, the exosomes are ready for downstream analysis. The exosome re-suspension is recommended to be stored at $-80 \text{ }^{\circ}\text{C}$ immediately.

Notice

This kit is for research use only, not for clinical diagnostic purpose.

$1 \times$ sterile PBS is not supplied and should be prepared by user.

We recommend that exosomes used for electron microscopy, NTA analysis and proteomics studies should be filtered by 0.22 μm filtration.

More detail information is on the official website: www.wayenbio.com

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