Wayen Exosome Isolation Kit

User Manual

Cat# EIQ3-02001 (Plasma) Version 2017-01

Description (For research only)

EIQ3-02001 kit is used to isolate / purify exosomes between 30 and 200 nm diameter from plasma. By adding appropriate amount of reagents to plasma sample, and incubating the mixture within a short period, the exosomes can be collected after sample centrifugation.

Advantages

- ✓ Quantity: Higher yield (versus other kits or methods)
- ✓ Quality: Pure exosome (less plasma high-abundant protein)
- ✓ Quick: Faster (< 2 hours)

Contents

EIQ3-02001 kit contains Reagent A, Reagent B and Reagent C.

Storage

The kits are shipped at 4 °C and should be stored properly after received. Properly stored kits are stable for 1 year from the date received.

Components	Storage	Amount
Reagent A	4 °C	7.5 mL
Reagent B	4 °C	7.5 mL
Reagent C	-20 °C	600 µL

Experiment Protocol of EIQ3-02001 (Plasma)

1. Prepare Sample

- 1.1 Take the plasma sample from storage and keep it on ice. If starting with frozen sample, thaw the sample completely in a 25 °C water bath and then place it on ice.
- 1.2 Take out Reagent C and thaw it completely on ice.

1.3 Add 4 μL Reagent C into 200 μL plasma and mix them well by vortexing or pipetting up and down until obtain a homogenous mixture.

※ **Note: Plasma : Reagent C = 50:1 (volume)**

1.4 Incubate the mixture at 37 $^{\circ}\text{C}$ for 15 minutes.

1.5 After incubation, the samples turn into jellylike status. Tap the tubes firmly to change them into liquid status and then centrifuge at $10,000 \times g$ for 10 minutes at room temperature.

1.6 Transfer the supernatant to a new 1.5 mL tube and then place it on ice.

2. Isolate Exosomes (Balance the Reagent A and Reagent B to room temperature before use and the starting volume of plasma is recommended to be 200 μL . The example below is shown with 200 μL plasma)

2.1 Take out 200 μL pre-treated plasma. Add 50 μL Reagent A into it and mix well by pipetting up and down or vortexing until obtain a homogenous mixture.

※ **Note: Plasma : Reagent A = 4 : 1 (volume)**

2.2 Incubate the mixture at 4 $^{\circ}\text{C}$ for 30 minutes.

※ **Note: Do not rotate or shake the tube during the incubation period**

2.3 After incubation, centrifuge the mixture at $3000 \times g$ for 10 minutes at room temperature. Remove the supernatant and the pellet is at the bottom of the tube.

2.4 Centrifuge the sample again within a short moment and remove the residual supernatant.

2.5 Resuspend the pellet completely by adding 200 μL 1 \times sterile PBS. Mix well by pipetting up and down or vortexing until obtain a homogenous mixture.

※ **Note: Volume of re-suspension is equal to the starting plasma volume at this procedure**

2.6 Add 50 μL Reagent B to the re-suspension and mix well by pipetting up and down or vortexing until obtain a homogenous solution.

※ **Note: The re-suspension : Reagent B = 4:1 (volume)**

2.7 Incubate the mixture at 4 $^{\circ}\text{C}$ for 30 minutes.

2.8 After incubation, centrifuge the mixture at $3000 \times g$ for 10 minutes at room temperature and remove the supernatant.

2.9 Centrifuge the sample again within a short moment and remove all residual supernatant.

2.10 Resuspend the exosome pellet completely in 50-120 μL 1 \times sterile PBS and mix well to obtain a homogenous mixture. Once the pellet is re-suspended, the exosomes are ready for downstream analysis. The exosome re-suspension is recommended to be stored at -80°C immediately.

Notice

This kit is for research use only, not for clinical diagnostic purpose.

1 \times sterile PBS is not supplied and should be prepared by user.

We recommend that exosomes used for electron microscopy, NTA analysis and proteomics studies should be filtered by 0.22 μm filtration.

More detail information is on the official website: www.wayenbio.com

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