Wayen Exosome Isolation Kit

User Manual

Cat# EIQ3-04001 (Cell Culture Media) Version 2017-01

Description (For research only)

EIQ3-04001 kit is used to isolate exosomes between 30 and 200 nm diameter from cell culture media. By adding appropriate amount of reagent to the cell culture media, and incubating the mixture over 12 h, the exosomes can be collected with centrifugation.

Advantages

- ✓ Quantity: Higher yield
- ✓ Easy: No ultracentrifugation

Contents

EIQ3-04001 kit contains one reagent.

Storage

The kits are shipped at room temperature and should be properly stored after receipt. Properly stored kits are stable for 1 year from the date received.

Components	Storage	Volume
Reagent	4 °C	50 mL

Experiment Protocol of EIQ3-04001 (Cell Culture Media)

1. Prepare Sample

1.1 Take the cell culture media from storage and keep it on ice. If starting with frozen sample, thaw the sample completely in a 25 °C water bath and then place it on ice.

1.2 Centrifuge the cell culture media at $3,000 \times g$ for 15 minutes at $4\text{ }^{\circ}\text{C}$.

1.3 Transfer the supernatant to a new tube and place it on ice until ready to perform the isolation.

2. Isolate Exosomes (Balance the Reagent to room temperature before use and the starting volume of cell culture media is recommended to be 10 mL. The protocol below is shown with 10 mL cell culture media.)

2.1. Take out 10 mL pre-treated cell culture media, add 5 mL Reagent into it and invert the mixture up and down until obtain a homogenous mixture.

※**Note: Cell culture media : Reagent =2:1 (volume)**

2.2 Incubate the mixture at $4\text{ }^{\circ}\text{C}$ overnight (12-16 h).

2.3 After incubation, centrifuge at $3,000 \times g$ for 60 minutes at $4\text{ }^{\circ}\text{C}$.

2.4 Take out 1 mL supernatant into a 1.5 mL tube firstly and remove the residual supernatant.

2.5 Resuspend the pellet completely with 1 mL supernatant above. Transfer the mixture to a new 1.5 mL tube.

2.6 Centrifuge the re-suspension at $10,000 \times g$ for 10 minutes at $4\text{ }^{\circ}\text{C}$. Remove the supernatant without disturbing the precipitated pellet.

2.7 Resuspend the pellet with 50-200 μL $1 \times \text{PBS}$, and mix it well by vortexing or pipetting up and down until obtain a homogenous mixture.

2.8 Centrifuge the sample again at $10,000 \times g$ for 5 minutes at $4\text{ }^{\circ}\text{C}$. Transfer the supernatant to a new tube.

2.9 The supernatant contains exosomes. The exosomes can be used for downstream analysis immediately or aliquoted and stored at $-80\text{ }^{\circ}\text{C}$ till next experiment.

Notice

This kit is for research use only, not for clinical diagnostic purpose.

$1 \times \text{PBS}$ is not supplied and should be prepared by user.

We recommend that exosomes used for electron microscopy, NTA analysis and proteomics studies should be filtered by $0.22\text{ }\mu\text{M}$ filtration.

For more detail information, please visit our official website: www.wayenbio.com

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